



**Klausz Dental
Laboratories Ltd.**

klauszdl@istar.ca
www.klauszdentallab.com

"Working Harder And Smarter For Your Practice"

40 Kodiak Crescent, Units 11-12 Toronto
ON Canada M3J 3G5
Toll Free: 1-800-410-4392
Local: (416) 630-6262

Quality | Technology | Innovation

Olmos Orthotics Indications




T&S Certified Lab

Olmos Night Appliance

ON1- Anterior Deprogrammer

Indication

- ◆ Control Parafunctional activity
- ◆ Reduce Clenching and Grinding
- ◆ Acute or chronic bruxism




Olmos Day Appliance

OD1- Acrylic Day

Indication-


- ◆ Indicated for chronic disc-displacement or skeletal and muscular asymmetry.
- ◆ Functions as a positioning stent to correct vertical dimension, rotation, cant and protrusive changes.
- ◆ Lingual Anatomical ridge guides disclusion and protects contralateral joint in working movements.



ON2- Night Positioner

Indication


- ◆ Prevent disc dislocation (locking) in the supine
- ◆ Protect inflamed or injured tissues



OD2- Compact Day

Indication- Same as OD1 with wire lingual connector for patients with:

- ◆ Narrow arches
- ◆ Anterior crowding
- ◆ Require minimal lingual interference



ON3- Night Positioner

Indication-

- ◆ Combines ON1 and ON2
- ◆ Reduces inflammation (capsulitis)
- ◆ Prevent locking in supine position
- ◆ Improve airway obstruction
- ◆ Reduces clenching or grinding



OD3- PMT Thermoform Day

Indication- Same as OD1 with thermoform design for patient with:

- ◆ Narrow arches
- ◆ Inter-occlusal distances less than 1.5mm
- ◆ Patients with allergies to metals
- ◆ Not recommended for patients with severe crowding or undercuts

