



**Klausz Dental
Laboratories Ltd.**

Working **harder and smarter** for your practice!

www.klauszdentallab.com • klauszdl@uniserve.com • Tel: 416-630-6262 • 1-800-410-4392
40 Kodiak Crescent • Units 11-12 • Toronto, ON • Canada M3J 3G5

An O.I.C.S. Specialty Laboratory
Occlusion - Implants - Cosmetics - Sleep

Bite Repositioning Appliances
**The Perfect Accompaniment to Your
Panthera D-SAD Sleep Appliance**

We are proud to announce that we are a leading provider of Panthera D-SAD Appliances to dentists in Ontario



After wearing a Mandibular Advancement Device (MAD) over a period of time, there is an inevitable undesirable bite change with the potential of a posterior open bite developing, as the mandible slowly migrates forward. To mitigate this problem, it is advisable to wear a Bite Repositioning Appliance immediately after removing the MAD. The patient should bite on this until he feels a normal, solid contact on the back teeth. The length of time needed for this varies from patient to patient. This could be done while showering or commuting.

These Bite Repositioning Appliances are made from a thin thermoform material that fits over the lower teeth, similar to a sports mouth guard. The posterior occlusion of the maxillary opposing dentition is imprinted into the appliance and gently guides and holds the mandible in place as the muscles reposition the mandible into the normal occlusion.




Bite Repositioning Appliances from Klausz Dental Laboratories are available in a variety of colours and are designed to be as individual as your patients that wear them.




Clear




Purple




Orange




Green



Red



Yellow



Blue



Black

Always looking for ways to improve the quality of our work is just another way we are, "Working Harder and Smarter for your Practice!"