

PANTHERA SLEEP ADJUSTMENT FOR D-SAD RECOMMENDED TOOLS AND COMPOUND

ADJUSTMENT BURRS

BRASSELER USA
H79DF.11.040 HP DF
Cutter Carbide




Cutting Length (mm): 4.2
Head Size (1/10 mm): 023
Package Label: H73EF.11.023
Carbide Cutter
Shank: 11-Straight Handpiece

BRASSELER USA
H79DF.11.040 HP DF
Cutter Carbide



Cutting Length (mm): 13.0
Head Size (1/10 mm): 040
Package Label: H79DF.11.040
Carbide Cutter 1P
Shank: 11-Straight Handpiece
Shape: Round-End Taper


BRASSELER USA
H79GE.11.050 HP GE
Cutter Carbide



Cutting Length (mm): 13.0
Head Size (1/10 mm): 050
Package Label: H79GE.11.050
Carbide Cutter 1P
Shank: 11-Straight Handpiece
Shape: Round-End Taper

ELASTIC GROOVES

RIO GRANDE USA LYNX™
Crown Saw Burs
#348513



Brand: LYNX™
Head shape: Crown saw
Shank size: 3/32"

RIO GRANDE USA LYNX™
High-Speed Steel Twist Drills
#349206



Brand: LYNX™
Material: High-speed steel
Head size: 0.6mm
Shank size: 3/32"

FINISHING & POLISHING

RIO GRANDE USA
Fabulustre Cut
and Polish Compound
#331015



RIO GRANDE USA
Miniature Felt Point,
#330172



Material: FeltHead
Shape: Point
Dimensions: 3/4"H x 1/4" dia.
Shank size: 3/32"
Density: Hard